

# Index

## Numerics

2% bleach solution, definition of 6-5

## A

accuracy of dispense

*See* dispense accuracy

acetone in washes, special pump-head tubing for 6-6

ActiveX interface for use with ControlMate device-control software 1-7, 4-72

air aspirating into syringes (troubleshooting topic) 8-7

air bubble, restoring in syringe 7-26

air bubbles in NanoFill tubing, getting rid of E-56–E-58

air gap volume, maximum, resetting 4-13

air gaps

recommended volumes for 7-10

when to use 7-6

Air parameter

*See* ASPIRATE screen parameters

air-supply gauge setting E-10

applications

cell culture 7-16–7-17

creating membrane arrays and dot blots 7-18

external, executing 4-50–4-51

for Hydra II 384 base unit 1-6

for Hydra II 96 base unit 1-6

for Hydra II 96 base unit w/ R/B syringes 1-6

selection of correct syringe size for 5-1

using Hydra II with vacuum manifold and manual positioners 7-18

who to contact with questions about xvii

Aspirate air gap volume, resetting maximum 4-13

Aspirate command 4-55–4-56

aspirate depth 4-28, 4-59

Aspirate operation 3-9

Aspirate overstroke 4-20, 4-55

*See also* priming syringes

ASPIRATE screen parameters 3-5

Aspirate spitback volume 4-20, 4-55

Aspirate with Overstroke/Dispense Specific mode

priming syringes in 7-9

programming and running a file in 7-3–7-6

transferring microliter volumes with 7-2

when to use 7-2

Aspirate with Pre Air Gap/Dispense All mode

priming syringes in 7-9–7-10

programming and running a file in 7-7–7-8

using with minimal sample volumes 7-6

when to use 7-6

Aspr button 3-7

automated operation of microdispenser

description of 4-1–4-73

flowchart of steps for 4-3

Automatic Syringe Wash Module

illustration of 1-4

needle guide splashed during AutoWash operation (troubleshooting topic) 8-8

overview of 1-9

replacing pump-head tubing for 6-31–6-35

setting up 2-7–2-15, E-10–E-15

wash reservoir, using in AutoWash procedures 4-64

AutoWash procedures

advanced techniques for 6-10

definition of 6-5

guidelines for 6-6

sequence, description of 4-64

solutions verified for use with pump Wash Module 6-6

Wash operation 3-4

axes, homing 4-60

## B

backlash

*See* spitback volume

bleach

solution, definition of 6-5

using in wash procedures for nucleic acid samples 6-8

using with microdispenser xviii, 6-8

bottle assembly

*See* Pressure Reservoir

buttons on front of microdispenser 3-7–3-8

## C

C.V.  
*See* coefficient of variance  
 calibrating  
   labware  
     *See* labware  
   stage  
     *See* stage  
 calibration table, creating or editing E-62–E-63  
 Call Procedure command 4-47  
 cautions  
   *See* warnings and cautions  
 cell culture applications, using Hydra II in  
   7-16–7-17  
 Check vessel holder sensor option 4-59, 4-69  
 chemicals compatible with Hydra II 7-25  
 chemicals routinely used with Hydra II 7-26  
 cleaning wires, using to clean needles 6-16–  
   6-18  
 coefficient of variance (C.V.), precision of dis-  
   pense volume across syringe array 5-3  
 comb plate difficult to insert (troubleshooting  
   topic) 8-7  
 command tree 4-42  
 commands  
   *See* communications protocol; ControlMate  
   software; *individual command names*;  
   file parameters; functions; *and* Hydra  
   II-Plus-One System, using software  
 Commands tab, viewing and setting default val-  
   ues on 4-15  
 communications ports, setting up 4-14, 4-25,  
   E-24  
 communications protocol D-2–D-18  
 compressed-air supply, connecting E-9–E-10  
 configuration, of software, checking 4-12  
 connecting the microdispenser to a computer  
   2-16  
 Connection tab, viewing to check hardware con-  
   nections 4-14  
 contamination between sample preparations,  
   preventing 7-13–7-16  
 controller unit, Hydra II-Plus-One E-9–E-10,  
   E-15  
 ControlMate OLE software component 1-7,  
   4-72

ControlMate software  
   *See also* labware; sequence files; *and* stage  
   actions within 4-11  
   ActiveX interface optional component for  
     4-72  
   child dialogs in 4-11  
   command groups in 4-10  
   command reference for Hydra II 4-39–4-68  
   command-parameter values, default, setting  
     4-15  
   commands, list of 4-39  
   help system in 4-6  
   installing 4-4  
   online help 4-2  
   overview of 4-4–4-11  
   system requirements for A-1  
 corrosive agents  
   *and* compatibility with Hydra II 7-26  
   protecting microdispenser from xviii  
 CPG beads, using with Hydra II 96 with R/B sy-  
   ringes 7-24  
 Cycles/Sec, setting in Nano Dispense function  
   E-39

## D

default command-parameter values, setting in  
   ControlMate 4-15  
 Device setting, checking 4-13  
 Disp button 3-8, 3-9  
 dispense accuracy  
   measuring with gravimetric test 5-4  
   measuring with standard curve 5-5–5-6  
 dispense accuracy and precision, for Hydra II-  
   Plus-One, validating E-59–E-61  
 Dispense command 4-57  
 dispense height  
   for very small dispense volumes 7-12  
   setting 3-10, 4-28, 4-59, 7-10–7-13  
   when dispensing into wet plates 7-10–7-11  
   when dispensing onto dry plates 7-11–7-13  
 Dispense operation 3-9  
 dispense precision  
   according to syringe size 5-1–5-2  
   evaluating 5-6–5-8  
 dispense problems with small volumes (trouble-  
   shooting topic) 8-5  
 DISPENSE screen parameters 3-5  
 dispense volumes, and pulse-width settings, for  
   Hydra II-Plus-One E-59–E-61

- dispensing into wet plates, dispense height for 7-10-7-11
  - dispensing onto dry plates
    - dispense height for 7-11-7-13
    - dispensing very small volumes 7-12
    - need to prevent "creep" at needle tip 7-26
    - uneven dispense (troubleshooting topic) 8-6
  - dispensing onto sidewalls of CyclePlate wells to prevent cross-contamination 7-14, 7-15
  - display, LCD, on microdispenser
    - displaying custom message on, from running sequence file 4-66
    - setting file parameters on 3-5-3-16
  - DMSO-based samples
    - preventing carry-over contamination in 7-14
    - special pump-head tubing for 6-6
    - wash procedures for 6-7-6-8
  - dot blots, creating 7-18
  - Down button 3-7, 3-9
  - DuraFlex needles 1-8
  - dwelt time 4-55, 4-57
- E**
- Edit button 3-8
  - Edit Lock 4-11
  - electrical interference xviii
  - electrical specifications for Hydra II system A-2-A-3
  - Emergency Shutoff Switch xv
    - installing 2-6, E-16
    - stopping the system with xv, 3-18, 4-69
    - testing 2-7, E-17
  - Empty button 3-7, 3-9
  - Empty operation 3-9
  - EMPTY screen parameters 3-6
  - EPROM file parameters
    - See* file parameters
  - error, in sequence-file run 4-14
  - Execute External Applications command 4-50-4-51
  - External Procedure File command 4-53
- F**
- fiber-optic sensors
    - checking xvi, 6-11
    - lenses, cleaning 6-11
  - file execution
    - See* sequence files, running
  - file parameters
    - setting (manual operation) 3-12-3-16
    - table of 3-5-3-7
  - file validation
    - See* sequence files, validating
  - files
    - See* file parameters *and* sequence files
  - finger slots on microdispenser sides 2-2
  - firmware version 3-13
  - fluid supply line, installing E-22
  - Fluidics Module
    - See* NanoFill Fluidics Module
  - functions
    - inserting in protocol
      - See* protocols, creating and editing
    - setting Cycles/Sec in Nano Dispense E-39
    - setting Tip Touch in Hydra Trial Dispense or Hydra Dispense E-38
    - setting Volume in Hydra Aspirate or Hydra Dispense E-37
  - fuses, checking and replacing 6-26-6-30
- G**
- glycerol in samples, wash procedures for 6-7-6-8
  - gravimetric test, using to measure dispense accuracy 5-4-5-5
  - Group Commands command 4-49
- H**
- hands, protecting xvi
  - hardware specifications for Hydra II systems A-2-A-3
  - hardware, installing 2-1-2-16, E-7-E-23
  - helium supply, for Hydra II-Plus-One E-7, E-22
  - help
    - online, in ControlMate software 4-2, 4-6
    - with questions about applications for microdispenser xvii
    - with service and maintenance tasks 6-3
  - Home Axes command 4-60
  - Hydra II 96 base unit, description of 1-6
  - Hydra II 96 with R/B syringes
    - description of 1-6
    - using resin suspensions with 7-22-7-24
  - Hydra II 96, using with manually operated plate positioners 7-19-7-22
  - Hydra II doesn't operate although power is on (troubleshooting topic) 8-3

- Hydra II EPROM file parameters
    - See file parameters
  - Hydra II files
    - See file parameters *and* sequence files
  - Hydra II is "behaving" erratically (troubleshooting topic) 8-2
  - Hydra II microdispenser
    - See also Hydra II-Plus-One System; service and maintenance; *and* system applications for xvii, 5-1, 7-16–7-19
    - automated operation of 4-1–4-73
    - buttons, layout of 3-7
    - chemicals compatible with 7-25
    - chemicals routinely used with 7-26
    - communications protocol for D-2
    - connecting to computer 2-16
    - dispense accuracy
      - See dispense accuracy
    - dispense precision
      - See dispense precision
    - flowchart for automated operation of 4-3
    - flowchart for manual operation of 3-2
    - illustration of 1-2
    - installing 2-1–2-16
    - interface, flowchart for manual operation 3-3
    - manual operation of 3-1–3-21
    - models available 1-5
    - moving 2-2
    - overview of 1-1
    - Plus-One option
      - See Hydra II-Plus-One System
    - protecting from damage xv, xvii–xix
    - quick tips for use of 8-8
    - related products for 1-10, B-4–B-5
    - replacement of packaging materials for B-3
    - returning for service or repair xi
    - setup of 2-1–2-16
    - site preparation for 2-1
    - stopping xv
    - surge protection for xviii
    - syringes for 1-8
    - troubleshooting
      - See troubleshooting topics
    - unpacking 2-1–2-2
    - using viscous solutions with 7-24–7-25
  - Hydra II microdispenser (continued)
    - validating performance of
      - See dispense accuracy *and* dispense precision
    - Wash Module
      - See Automatic Syringe Wash Module
    - X/Y Plate Stage
      - See stage
  - Hydra II operations, executing with microdispenser buttons 3-9–3-10
  - Hydra II Wash operation, flowchart for 3-4
  - Hydra II with Automatic Syringe Wash Module
    - See also Automatic Syringe Wash Module
    - description of 1-7
    - illustration of 1-4
  - Hydra II with X/Y Plate Stage
    - See also stage
    - description of 1-7
    - illustration of 1-4
  - Hydra II won't turn on (troubleshooting topic) 8-2
  - Hydra II-Plus-One System
    - See also NanoFill dispenser; NanoFill Fluidics Module; *and* protocols
    - components of E-2–E-3
    - description of 1-7
    - overview of E-1–E-5
    - performance, validating E-59–E-61
    - priming E-56–E-58
    - protocols E-34–E-56
    - roadmap of system tasks E-4
    - setting up E-7–E-23
    - setup, quickview of E-4–E-6
    - stopping E-51
    - troubleshooting E-56
    - turning on, order of steps for E-23
    - using software E-23–E-56
  - Hydra-96 Plate Positioner for CyclePlate-384 and -192, using 7-22
  - Hydra-96 Plate Positioner for Standard 384 Well Plates, using 7-21
- I**
- Instrument Return Report Form, when to use xi
  - IRR Form
    - See Instrument Return Report Form
- L**
- labware
    - See also dispense height

labware (continued)

creating, editing, and deleting definitions for 4-26–4-34, E-28–E-33

setting up Tip Touch on Hydra II-Plus-One E-30

types available for different systems 4-58

types of and required parameters for 4-30

using manually operated plate positioners with 7-19–7-22

well-center-to-well-center spacing in 7-20, 7-21

lash, accounting for by priming syringes 7-9

Limits tab, viewing and setting default values on 4-15

linear regression curve (LRC), using to validate system performance E-59–E-61

liquid-level sensor for Wash Reservoir 2-14

logic power supply board, changing fuses on 6-28

loops in protocols

*See* protocols, setting iterations for, on Hydra II-Plus-One

LRC

*See* linear regression curve

**M**

Main Application screen 4-5–4-10

Main Sequence Editor 4-39–4-45

manual operation of microdispenser

description of 3-1–3-21

flowchart for 3-2

manual wash procedures

guidelines for 6-7

on system with Wash Module 3-16, 4-64

on system without Wash Module 4-64

membrane arrays, creating 7-17–7-18

Menu Bar 4-6

Message Display command 4-66

microdispenser interface, manual operation, flowchart for 3-3

microliter volumes, transferring 7-2

microplates

*See* labware

mineral oil or glycerol in samples, wash procedures for 6-7–6-8

Mix command 4-65

Mix operation, using Wash operation as 7-19

motor not working (troubleshooting topic) 8-3

motor power supply board, changing fuses on 6-28

Motors tab, viewing settings on 4-17

Move to Position command 4-58–4-59

moving the system 2-2

**N**

Nano Dispense Custom protocols, on Hydra II-Plus-One E-47–E-52

NanoFill dispenser

installing E-19–E-20

leaking, how to fix E-58

manual run E-53–E-56

priming E-50

recalibrating if necessary E-62–E-63

validating performance of E-59–E-61

washing E-64

NanoFill Fluidics Module

connecting tubing and hookups for E-21–E-22

fluid supply line, installing E-22

setting up E-18–E-19

wash reservoir for E-20

needle guide, correct position for xviii

needles

cleaning 6-14–6-25

cleaning clogged needles with cleaning wires 6-16–6-18

design of 1-8

not aligning with microplate (troubleshooting topic) 8-4

preventing "creep" at tips of 7-26

protecting xviii

solutions wicking up, how to prevent 7-26

types available 1-8

nest sensors

turning on or off 4-59

Notepad command 4-48

nucleic acid samples, wash procedures for 6-7–6-8

**O**

offset values in labware definitions, applying 4-28

online help, in ControlMate 4-2

operations screens, buttons on 3-8

operator, protecting xvi–xvii

organic solvents, wash procedure for 6-7

overstroke

*See* Aspirate overstroke

- oxidizing agents
  - and compatibility with Hydra II 7-26
  - protecting microdispenser from xviii
- P**
- P1Time parameter
  - See Pumps screen parameters
- P2Time parameter
  - See Pumps screen parameters
- P3Time parameter
  - See Pumps screen parameters
- packaging materials replacement B-3
- packing list xiii
- Page button 3-7
- passwords for sequence files 4-11
- Pause command 4-61
- PCR reactions, preventing carry-over contamination in 7-13
- performance of Hydra II microdispenser
  - See dispense accuracy *and* dispense precision
- physical specifications for Hydra II system
  - A-2-A-3
- Pipettor setting, checking 4-13
- plate positioners, manually operated, using with Hydra II 96 1-10, 7-19-7-22
- plates
  - See labware
- Plates tab, using to create, edit, and delete labware definitions 4-27-4-34
- Plus-One system
  - See Hydra II-Plus-One System
- pneumatic Z controller
  - See controller unit
- pooling multiple samples 7-6
- ports, communication, setting up E-24
- Positions tab, using to calibrate stage 4-20
- power cord for microdispenser, installing 2-6, E-16
- power module
  - fuses, checking and replacing 6-26-6-28
  - location of 2-4
- power switch, stopping the system with xv
- precision of dispense
  - See dispense precision
- Pressure Reservoir, setting up E-20-E-21
- preventing carry-over contamination
  - for DMSO-based samples 7-14
  - in sequencing reactions 7-14
  - when preparing PCR reactions 7-13
- preventing contamination between sample preparations 7-13-7-16
- preventing cross-contamination by dispensing onto sidewalls of CyclePlate wells 7-14
- preventing solutions from wicking up the sides of needles 7-26
- Prime parameter
  - See ASPIRATE screen parameters
- priming NanoFill dispenser E-50
- priming syringes 4-20, 7-8-7-10
- priming the system
  - See priming syringes
- printed circuit board
  - and grounding xviii
  - electrical interference affecting xviii
- Procedure command 4-52
- Procedure Library command 4-46
- procedures (ControlMate)
  - calling 4-47
  - external 4-53
- procedures (task lists)
  - calibrating the stage 4-19-4-24
  - checking and replacing motor power supply and logic power supply fuses 6-28-6-30
  - checking and replacing power-module fuses 6-26-6-28
  - cleaning needles with cleaning wires 6-16-6-18
  - cleaning needles with wash and abrasive material 6-15-6-16
  - creating a sequence file 4-36
  - creating, editing, and deleting labware 4-30-4-34
  - editing a sequence file 4-37-4-38
  - replacing pump-head tubing 6-31-6-35
  - replacing syringes 6-19-6-25
  - replacing Teflon tips on plungers 6-13-6-14
  - running a procedure 3-18-3-21
  - running a sequence file 4-69-4-72
  - setting file parameters 3-12-3-16
  - setting tray table/stage height (manual operation) 3-10
  - setting up manual wash on system with Wash Module 3-16
  - testing fiber-optic sensors 6-11
  - washing syringes at the end of the day 6-9-6-10
  - washing syringes between samples 6-7-6-9

products for Hydra II microdispenser 1-10  
 program, starting E-23  
 programming  
   *See* file parameters *and* sequence files  
 programs  
   *See* file parameters *and* sequence files  
 protecting the Hydra II microdispenser from  
 damage xvii  
 protecting the operator xvi–xvii  
 protein buildup in syringes, removing 6-9  
 protocols  
   and using Tip Touch, on Hydra II-Plus-One  
   E-38  
   creating and editing, on Hydra II-Plus-One  
   System E-34–E-49  
   Hydra II-Plus-One E-34–E-56  
   Nano Dispense Custom, creating, on Hydra  
   II-Plus-One E-47–E-49  
   running, on Hydra II-Plus-One E-50–E-52  
   setting iterations for, on Hydra II-Plus-One  
   E-44  
   stopping run of, on Hydra II-Plus-One E-51  
 pulse-width settings, converting to dispense vol-  
 umes E-59, E-59–E-61  
 Pump Control command 4-67  
 pump-head tubing on Wash Module, replacing  
 6-31–6-35  
 PUMPS screen parameters 3-7  
 pumps, controlling without use of wash routine  
 4-67  
 PWUS  
   *See* pulse-width settings

## R

Ready Mode screen File parameter 3-5  
 Ready Mode screen, buttons on 3-7  
 refurbishing syringes B-2  
 reservoirs  
   *See* labware  
 resin suspensions  
   examples of types used with Hydra II 7-23  
   using with Hydra II 96 with R/B syringes  
   7-22–7-24  
 returning microdispenser or microdispenser  
 parts for service or repair  
   instructions xi  
   packaging materials for B-3  
 roadmap of system tasks for Hydra II-Plus-One  
 E-4  
 run error, in sequence files 4-14

Run Locks 4-11

## S

safety  
   and position of needles xvii, xviii  
   and using sodium azide xvi, 6-11  
   checking fiber-optic sensors xvi, 6-11  
   keeping syringes from damage xvii  
   reporting uses of microdispenser on Instru-  
   ment Return Report Form xi  
   rinsing pump tubing before replacing 6-32  
   unplug Hydra II before checking motor  
   power supply or logic power supply fus-  
   es 6-28  
   using dispenser for intended applications  
   xvii  
   warnings and cautions xv–xvii  
 sample preparations, preventing contamination  
   between 7-13–7-16  
 sample volumes, very small, handling 7-6  
 samples, integrity of xix, 4-70  
 sensors  
   *See also* fiber-optic sensors  
   liquid-level sensor for Wash Reservoir 2-14  
   nest, turning on or off 4-59  
 Sephadex beads, using with Hydra II 96 with R/  
 B syringes 7-24  
 Sequence File Editor  
   *See* Main Sequence Editor  
 sequence files  
   creating and editing 4-35–4-68  
   executing  
     *See* sequence files, running  
   locking and unlocking 4-11  
   passwords for 4-11  
   rules for 4-38  
   run error 4-14  
   running 4-69–4-72  
   sample C-1–C-33  
   saving 4-36  
   stopping or pausing run of 4-72  
   testing 4-69  
   validating 4-36  
 sequencing reactions, preventing carry-over  
 contamination in 7-14  
 serial connection, testing 4-14  
 service and maintenance  
   checking and replacing fuses 6-26–6-30  
   checking fiber-optic sensors 6-11  
   checking Teflon tips 6-12

- service and maintenance (continued)
  - cleaning clogged needles with cleaning wires 6-16–6-18
  - cleaning needles 6-14–6-18
  - prolonging life of syringes and Teflon tips 6-9
  - replacement guidelines for Teflon tips 6-12
  - replacing pump-head tubing on Wash Module pump 6-31–6-35
  - replacing syringes 6-19–6-25
  - replacing Teflon tips 6-12–6-14
  - returning microdispenser or microdispenser parts for service or repair xi, 6-3
  - storing syringes for more than a week 6-11
  - table of service and maintenance tasks 6-2
  - washing NanoFill dispenser E-64
  - washing syringes 6-3–6-10
  - when to call Customer Service 6-1
- setting tray table/stage height 3-10, 4-28, 4-59
- setting up the system 2-1–2-15
- shipping plates, removing from stage 2-3, E-8
- shutting down the system xv, 3-18, 4-69
- single dispenser
  - See* NanoFill dispenser
- small volumes, transferring 7-2
- sodium azide, safe use of 6-11
- software
  - See* ControlMate software *and* Hydra II-Plus-One System, using software
- solid supports for chemistry, using with Hydra II 96 with R/B syringes 7-23
- SPA beads, using with Hydra II 96 with R/B syringes 7-23
- specifications for system A-1–A-3
- Speed Control command 4-68
- Speed parameter
  - See* ASPIRATE, DISPENSE, EMPTY, *and* WASH screen parameters
- spitback volume, resetting 4-13
- stage
  - attaching to microdispenser 2-2, E-8
  - calibrating 4-18–4-24, E-26–E-27
  - description of 1-7
  - height
    - See* setting tray table/stage height
  - homing before sequence-run file run 4-71
  - illustration of 1-5
  - protecting xviii, 2-3, E-9
  - recalibrating if moved 2-2
- stage (continued)
  - removing shipping plates from 2-3, E-8
- stainless steel needles 1-8
- standard curve, using to measure dispense accuracy 5-5–5-6
- Step 1 parameter
  - See* Pumps screen parameters
- Step 2 parameter
  - See* Pumps screen parameters
- sterile environment for transfer of cells and cell culture media 7-16
- STOP button 3-8
- stopping the system
  - with the Emergency Shutoff Switch xv, 3-18, 4-69, E-51
  - with the power switch xv
- surge protection for Hydra II xviii
- syringe array
  - changing from one size to another 2-16
  - mapping performance of 5-6
  - measuring precision and linearity of dispense volumes for 5-3–5-4
- syringes
  - See also* NanoFill dispenser
  - air aspirating into (troubleshooting topic) 8-7
  - avoiding damage to xvii
  - checking for damage to (troubleshooting topic) 8-6
  - clogged, how to clean 6-16–6-18
  - design and construction of 1-8
  - dispense precision according to syringe size 5-1
  - dispense problems with (troubleshooting topic) 8-5–8-6
  - emptying before running a sequence file 4-71
  - evaluating & mapping performance of 5-6
  - priming 4-20, 7-8–7-10
  - prolonging life of 6-9
  - protecting from damage xvii
  - refurbishing B-2
  - replacement parts for B-1–B-2
  - replacing 6-19–6-25
  - restoring air bubble in 7-26
  - sample wash procedures for 6-6–6-10
  - selecting correct size for application 5-1
  - storing for more than a week 6-11
  - types of 1-8
  - washing 6-3–6-10

system

*See also* Hydra II microdispenser *and* Hydra II-Plus-One System  
 carrying handles on 2-2  
 changing voltage setting for 2-4  
 components 1-1-1-7  
 connections  
   checking 4-12  
   illustration of 2-6, E-16  
   testing 4-14  
 electrical specifications for A-2-A-3  
 finger slots on 2-2  
 hardware specifications for A-2-A-3  
 moving 2-2  
 physical specifications for A-2-A-3  
 priming  
   *See* syringes, priming  
 setup 2-1-2-15  
 shutdown xv, 3-18, 4-69  
 software requirements for A-1  
 stopping xv, 3-18, 4-69  
 unpacking 2-1-2-2  
 system shutdown xv, 3-18, E-51

**T**

Table parameter

*See* ASPIRATE, DISPENSE, EMPTY, *and* WASH screen parameters

Teflon tips

checking 6-12  
 not replacing if reddish-brown precipitate appears in syringes 6-12  
 prolonging life of 6-9  
 replacement guidelines for 6-12  
 replacing 6-12-6-14  
 seal damaged (troubleshooting topic) 8-5

Tip Touch

calibrating labware for, on Hydra II-Plus-One E-30  
 in Aspirate command, in ControlMate 4-55  
 in Dispense command, in ControlMate 4-57  
 using in protocols, on Hydra II-Plus-One E-38

Tip wash vessel

*See* Wash Reservoir

tip/well offset

positioning needles in CyclePlate wells 7-15  
 setting in Move to Position command 4-59

tip/well offset (continued)

  setting offset values on Positions tab 4-22  
   values, applying 4-29  
 Toolbar icons 4-7  
 Toolbox commands 4-8  
 transferring individual samples of microliter volumes 7-6  
 transferring microliter volumes 7-2  
 transferring small volumes, problem with uneven dispense (troubleshooting topic) 8-5  
 tray table  
   *See also* troubleshooting topics  
   description of 1-1  
   guidelines for when to remove 6-17  
   height, setting 3-10, 4-28, 4-59  
   homing before sequence-file run 4-71  
   realigning 2-8  
   slotted, installing 2-8  
 tray table doesn't move (troubleshooting topic) 8-4  
 tray table doesn't move to set height (troubleshooting topic) 8-4  
 tray table/stage height, setting 3-10, 4-28, 4-59  
 troubleshooting topics  
   air is aspirating into syringes 8-7  
   checking for syringe damage 8-6  
   comb plate difficult to insert 8-7  
   damage to seal on Teflon tips 8-5  
   Hydra II doesn't operate although power is on 8-3  
   Hydra II is "behaving" erratically 8-2  
   Hydra II won't turn on 8-2  
   motor not working 8-3  
   needle guide splashed during AutoWash operation 8-8  
   needles don't line up properly in wells 8-4  
   problems with dispensing small volumes 8-5  
   syringes not dispensing properly 8-5  
   tray table doesn't move 8-4  
   uneven dispense when dispensing onto dry plates 8-6  
   uneven dispense when transferring small volumes 8-5  
 troubleshooting, on Hydra II-Plus-One System E-56  
 tubing and hookups for NanoFill Fluidics Module, connecting E-21-E-22  
 tubing, compressed-air supply on Hydra II-Plus-One, connecting E-9-E-10

## U

uneven dispense when dispensing onto dry plates (troubleshooting topic) 8-6  
 uneven dispense when transferring small volumes (troubleshooting topic) 8-5  
 unpacking the system 2-1-2-2, E-7  
 Up button 3-7, 3-9

## V

Valve Control command 4-54  
 vessels  
   *See* labware  
 viscous solutions, use of Hydra II microdispenser with 7-24-7-25  
 Vol parameter  
   *See* ASPIRATE, DISPENSE, *and* WASH screen parameters  
 voltage setting for system 2-4

## W

warnings and cautions  
   avoiding contact with sodium azide 6-11  
   avoiding damage to pump-head tubing 6-35  
   protecting the operator xvi-xvii  
   protecting the system from damage xvii-xix  
   rinsing pump-head tubing before replacing 6-32  
   stopping the system xv-xvi  
   Teflon tips, when NOT to replace 6-12  
   unplugging microdispenser before checking fuses for motor power supply or logic power supply 6-28  
   using corrosive or oxidizing agents with microdispenser 6-4  
   validating wash procedures 6-3  
 Wash button 3-10  
 wash cycles  
   definition of 6-5  
   in manual operation  
     *See* WASH screen parameters  
 Wash Module  
   *See* Automatic Syringe Wash Module  
 Wash Needles command 4-62-4-64  
 Wash operation  
   description of (manual operation) 3-10  
   flowchart for 3-4  
   using as Mix operation 7-19

Wash parameter  
   *See* WASH screen parameters  
 wash procedures  
   at end of day 6-9-6-10  
 AutoWash  
   advanced techniques for 6-10  
   definition of 6-5  
   between samples 6-7-6-9  
   for DMSO-based samples 6-7-6-8  
   for NanoFill dispenser E-64  
   for nucleic acid samples 6-7-6-8  
   for organic solvents 6-7-6-8  
   for samples with mineral oil or glycerol 6-7-6-8  
   guidelines for 6-4, 6-6, 6-7  
   manual wash on system with Wash Module 3-16, 4-64  
   manual, definition of 6-5  
   sample 6-6-6-10  
   solutions verified for use with pump Wash Module 6-6  
   terminology for 6-5  
 wash pumps, controlling without use of wash routine 4-67  
 Wash Reservoir  
   *See also* Automatic Syringe Wash Module use with Wash Needles command 4-64  
   wash reservoir for NanoFill Fluidics Module, setting up E-20  
 WASH screen parameters 3-6  
 wash system  
   *See* Automatic Syringe Wash Module  
 wash volume, definition of 6-5  
 well-center-to-well-center spacing in microplates 7-20  
 wires, for cleaning needles 6-16

## X

X/Y Plate Stage  
   *See* stage  
 X-axis, homing 4-60

## Y

Y-axis, homing 4-60

## Z

Z-axis, homing 4-60